

luna

new track lifts, new ideas, new feelings





luna

lifting with care

With the luna track lift Ergolet raises technical aids to a higher level

Using a balanced blend of new technology, new ideas and with focus on design we have developed a track lift system that will revolutionize technical aids globally.

Track lifts are the future, giving the carer better working conditions and the client a safe and dignified transfer.

Track lift systems will be mobile and become an even more viable alternative to traditional floor mounted patient lifts. That's why we have developed luna.

Three lifts in one

A new and patented system enabling use in



Vertically mounted
- the lift is slim and
allows the client
more headroom

Luna is easy to fit and remove thus allowing use in adjacent rooms or other locations.

The luna mounted vertically gives the client a greater sense of security and comfort



Mounted
horizontally - the lift
is flat allowing
optimal transfer
height

Luna can also be mounted horizontally thus giving maximum lifting height and better working conditions for example through doorways and room to room transfer



Fixed - mounted
without using tools

Luna is easily mounted directly onto the track, safely, simply and totally without using tools.

Mounted as a traditional track lift, the new easy to use spreader bar combined with simple and safe functions make client transfer optimal.



Better working conditions

Functionality and psychological well being have long been overlooked in client care. The working environment for the carer is also in need of improvement - we would like to introduce this improvement with luna.

Room to room transfer made easy



Using two lifts

Transfer from room to room is safe and simple using two lifting units and one spreader bar.



Or just one lift

Transfer from room to room easily made using just one lift unit.



Client confidence and security

Our mission is to develop products that give our clients confidence and address the natural suspicion that many people have against machines. We focus therefore on small details that make big differences.



Increased safety with the patented spreader bar



New spreader bar system with spring loaded hooks

The transfer situation requires great confidence between client and carer.

The simple spring loaded hook safeguards against accidental unhooking of the sling lifting straps whilst the click lock mechanism ensures greater security for the client



A hand holding a green and silver hand control device. The device has a small screen at the top displaying '130 cm' and '24 v 9 kg'. Below the screen are several circular buttons with icons: a four-way directional pad, a 'Menu' button, a button with a vertical double-headed arrow, a button with a red wavy line, a button with a key icon, and a large red button with a white triangle. The brand name 'ergolet' and the number '612' are printed at the bottom of the device. A coiled black cable is attached to the bottom. The background is a blurred image of a person's face.

Hand control of the future

Emergency functions built into the hand control

The luna hand control is built with the future in mind. A large display and simple to use functions ease the day to day use whilst giving back up technicians new possibilities for optimal maintenance. The hand control has three models with the one shown being the most advanced.

Different people


Different spreader bars



Mount a different spreader bar - push in, twist and lock using tools


2, 3 and 4 point spreader bars in different widths ensure dignified and safe transfers in all situations.

Luna is fitted with double security, the lift cannot be operated unless the spreader is correctly fitted.



New Mobile track systems


The track lift becomes a reusable technical aid



Patented X_Y system

Full room coverage can be fitted in every room. The system, comprising three telescopic tracks and four corner posts, can be installed in under half an hour by two persons. No tools are necessary for installation.

The mobile X - Y transfer system allows track lifts to become reusable technical aids for the first time.




Single track system

A single track system is in the same way set up easily and quickly.


Mobile systems give huge advantages in the working situation since track lifts can be used in even more transfer and lifting situations. Track lifts become mobile technical aids.

Fixed track systems



Patented X-Y system

Full room coverage with tracks mounted on the ceiling or wall. The rounded track profiles makes it easier for the system to blend into the room environment.



Single track system

The fixed single track system is the most economical using both straight and curved tracks

luna

Lifting with care



Data and specifications

Arguments for luna

Meets every transfer demand

A light lift unit that can easily be shared between several locations

Emergency functions on both the lift unit and hand control

Charging via the hand control in the charger or directly on the lift unit

Mount and dismount without using tools

Spreader bar with unique lifting hook system

Room to room transfer

Portable for use in other locations

Used with flexible track systems to meet every need

Technical specifications

Lifting capacity: 200kg

Weight: 6,7kg

Batteries: 36V, lead acid batteries

ergolet
lift the world

Manufacturer:

Ergolet A/S
Tårnborgvej 12C
4220 Korsør
Denmark

Tel: 58 37 37 20
Fax: 58 37 31 72
E-mail: ergolet@ergolet.dk

Distributed in Canada by:

J.W. Westman Inc.
Phone: (800) 387-8204
Fax: (905) 821-3168
email: jwwestman@bellnet.ca
www.jwwestman.com